



**EIGHTCHECK-3WP X-TRA™**  
HEMATOLOGY CONTROL

poch-100i

Parameters	Low Abnormal			Normal			High Abnormal		
	Mean	Limit	BCR	Mean	Limit	BCR	Mean	Limit	BCR
Quality Control Data Due on the 20 <sup>th</sup> of each month									
Expiration Date: 2010-08-11									
Lot Number: 01240710									
WBC x10 <sup>3</sup> /μL	2.8	0.4		6.8	0.7		17.6	1.8	
RBC x10 <sup>6</sup> /μL	2.36	0.21		4.48	0.27		5.24	0.26	
HGB g/dL	6.2	0.3		13.2	0.7		16.9	0.7	
HCT %	19.0	1.6		38.8	2.7		48.3	3.9	
MCV fL	80.5	12.8		86.6	10.6		92.2	11.4	
MCH pg	26.3	3.4		29.5	3.1		32.3	2.8	
MCHC g/dL	32.6	4.0		34.0	3.8		35.0	3.9	
PLT x10 <sup>3</sup> /μL	47	16		192	35		450	59	
Lymph (W-SCR) %	23.8	3.6		30.6	3.1		35.5	3.6	
Mixed (W-MCR) %	9.9	4.5		10.8	2.7		15.0	2.6	
Neut (W-LCR) %	66.3	5.3		58.6	3.5		49.5	4.0	
Lymph (W-SCC) x10 <sup>3</sup> /μL	0.7	0.2		2.1	0.5		6.2	0.7	
Mixed (W-MCC) x10 <sup>3</sup> /μL	0.3	0.2		0.7	0.2		2.6	0.5	
Neut (W-LCC) x10 <sup>3</sup> /μL	1.8	0.4		4.0	0.6		8.8	1.1	
W-SMV fL	67.2	16.8		66.2	14.6		66.0	13.2	
W-LMV fL	211.3	21.1		211.9	21.2		212.4	21.2	
RDW-SD fL	41.2	7.0		43.0	6.9		47.6	7.1	
RDW-CV %	14.4	2.3		14.2	2.4		14.7	2.6	
MPV fL	10.1	1.6		9.8	1.1		9.7	0.9	

Refer to package insert for Sysmex recommendations on setting QC targets and limits  
Distributed by  
Sysmex America, Inc.