



EIGHTCHECK-3WP X-TRA™
HEMATOLOGY CONTROL

pocH-100i

Low Abnormal	
Insight Due Date #1	10/24/2011
Expiration Date:	2011-12-28
Lot Number:	12630710

Normal	
Insight Due Date #2	11/26/2011
Expiration Date:	2011-12-28
Lot Number:	12630711

High Abnormal	
Insight Due Date #3	12/28/2011
Expiration Date:	2011-12-28
Lot Number:	12630712

Parameters	Mean	Limit	BCR
WBC x10 ³ /μL	2.8	0.4	
RBC x10 ⁶ /μL	2.36	0.21	
HGB g/dL	6.0	0.3	
HCT %	18.4	1.6	
MCV fL	78.0	12.7	
MCH pg	25.4	3.4	
MCHC g/dL	32.6	4.1	
PLT x10 ³ /μL	50	17	
Lymph (W-SCR) %	24.8	3.7	
Mixed (W-MCR) %	10.5	4.6	
Neut (W-LCR) %	64.7	5.2	
Lymph (W-SCC) x10 ³ /μL	0.7	0.2	
Mixed (W-MCC) x10 ³ /μL	0.3	0.2	
Neut (W-LCC) x10 ³ /μL	1.8	0.4	
W-SMV fL	65.3	16.6	
W-LMV fL	198.8	20.5	
RDW-SD fL	44.1	7.1	
RDW-CV %	16.0	2.4	
MPV fL	10.3	1.6	

Parameters	Mean	Limit	BCR
WBC x10 ³ /μL	7.0	0.7	
RBC x10 ⁶ /μL	4.26	0.26	
HGB g/dL	12.1	0.7	
HCT %	35.9	2.7	
MCV fL	84.3	10.6	
MCH pg	28.4	3.1	
MCHC g/dL	33.7	3.8	
PLT x10 ³ /μL	193	36	
Lymph (W-SCR) %	31.0	3.1	
Mixed (W-MCR) %	10.9	2.8	
Neut (W-LCR) %	58.1	3.5	
Lymph (W-SCC) x10 ³ /μL	2.2	0.5	
Mixed (W-MCC) x10 ³ /μL	0.8	0.2	
Neut (W-LCC) x10 ³ /μL	4.0	0.9	
W-SMV fL	65.1	14.5	
W-LMV fL	199.4	20.3	
RDW-SD fL	43.3	6.8	
RDW-CV %	14.4	2.3	
MPV fL	10.1	1.1	

Parameters	Mean	Limit	BCR
WBC x10 ³ /μL	17.7	1.8	
RBC x10 ⁶ /μL	5.28	0.26	
HGB g/dL	16.8	0.7	
HCT %	48.4	3.9	
MCV fL	91.7	11.3	
MCH pg	31.8	2.7	
MCHC g/dL	34.7	3.9	
PLT x10 ³ /μL	464	60	
Lymph (W-SCR) %	35.8	3.5	
Mixed (W-MCR) %	15.7	2.7	
Neut (W-LCR) %	48.5	3.9	
Lymph (W-SCC) x10 ³ /μL	6.3	0.7	
Mixed (W-MCC) x10 ³ /μL	2.8	0.5	
Neut (W-LCC) x10 ³ /μL	8.6	1.1	
W-SMV fL	64.0	13.1	
W-LMV fL	199.2	20.3	
RDW-SD fL	46.8	6.8	
RDW-CV %	14.4	2.5	
MPV fL	9.9	0.9	

Refer to package insert for Sysmex recommendations on setting QC targets and limits
Distributed by
Sysmex America, Inc.