



**EIGHTCHECK-3WP X-TRA™**  
HEMATOLOGY CONTROL

pocH-100i

Parameters	Low Abnormal			Normal			High Abnormal				
	Mean	Limit	BCR	Mean	Limit	BCR	Mean	Limit	BCR		
	Insight Due Date #1 1/15/2012			Insight Due Date #2 2/17/2012			Insight Due Date #3 3/21/2012				
	Expiration Date: 2012-03-21			Expiration Date: 2012-03-21			Expiration Date: 2012-03-21				
Lot Number: 13470710			Lot Number: 13470711			Lot Number: 13470712					
WBC x10 <sup>3</sup> /μL	3.3	0.4		6.6	0.7		17.6	1.8			
RBC x10 <sup>6</sup> /μL	2.36	0.21		4.39	0.26		5.19	0.26			
HGB g/dL	6.3	0.3		13.5	0.7		16.7	0.7			
HCT %	18.8	1.6		38.6	2.7		47.1	3.9			
MCV fL	79.7	12.7		87.9	10.6		90.8	11.3			
MCH pg	26.7	3.4		30.8	3.1		32.2	2.7			
MCHC g/dL	33.5	4.1		35.0	3.8		35.5	3.9			
PLT x10 <sup>3</sup> /μL	51	17		194	36		448	60			
Lymph (W-SCR) %	23.7	3.7		31.1	3.1		39.6	3.5			
Mixed (W-MCR) %	12.2	4.6		10.2	2.8		14.0	2.7			
Neut (W-LCR) %	64.1	5.2		58.7	3.5		46.4	3.9			
Lymph (W-SCC) x10 <sup>3</sup> /μL	0.8	0.2		2.1	0.5		7.0	0.7			
Mixed (W-MCC) x10 <sup>3</sup> /μL	0.4	0.2		0.7	0.2		2.5	0.5			
Neut (W-LCC) x10 <sup>3</sup> /μL	2.1	0.4		3.8	0.9		8.1	1.1			
W-SMV fL	67.6	16.6		65.8	14.5		65.6	13.1			
W-LMV fL	201.4	20.5		199.2	20.3		199.5	20.3			
RDW-SD fL	40.0	7.1		43.3	6.8		46.8	6.8			
RDW-CV %	14.1	2.4		13.8	2.3		14.5	2.5			
MPV fL	10.7	1.6		9.7	1.1		9.5	0.9			

Refer to package insert for Sysmex recommendations on setting QC targets and limits  
Distributed by  
Sysmex America, Inc.