































































EIGHTCHECK-3WP X-TRA™
HEMATOLOGY CONTROL

poch-100i

Parameters	Low Abnormal			Normal			High Abnormal		
	Mean	Limit	BCR	Mean	Limit	BCR	Mean	Limit	BCR
	Data Due Date #1: 04-Nov-08 Expiration Date: 2008-12-31 Lot Number: 82670710 			Data Due Date #2: 02-Dec-08 Expiration Date: 2008-12-31 Lot Number: 82670711 			Data Due Date #3: 05-Jan-09 Expiration Date: 2008-12-31 Lot Number: 82670712 		
WBC x103/mL	2.3	0.3		6.5	0.7		16.5	1.7	
RBC x106/mL	2.38	0.21		4.43	0.27		5.24	0.26	
HGB g/dL	6.3	0.3		13.3	0.7		16.7	0.7	
HCT %	18.8	1.6		38.4	2.7		47.3	3.8	
MCV fL	79.0	12.6		86.7	10.6		90.3	11.2	
MCH pg	26.5	3.4		30.0	3.1		31.9	2.8	
MCHC g/dL	33.5	4.2		34.6	3.8		35.3	3.9	
PLT x103/mL	50	17		194	35		469	61	
Lymph (W-SCR) %	23.8	3.6		30.1	3.0		35.1	3.5	
Mixed (W-MCR) %	6.8	3.1		10.7	2.7		14.5	2.5	
Neut (W-LCR) %	69.4	5.6		59.2	3.6		50.4	4.0	
Lymph (W-SCC) x103/mL	0.5	0.2		2.0	0.5		5.8	0.7	
Mixed (W-MCC) x103/mL	0.2	0.1		0.7	0.2		2.4	0.5	
Neut (W-LCC) x103/mL	1.6	0.3		3.8	1.6		8.3	1.0	
W-SMV fL	66.4	16.6		66.2	14.6		66.2	13.2	
W-LMV fL	209.6	21.0		209.1	20.9		209.1	20.9	
RDW-SD fL	45.4	7.7		45.6	7.3		45.2	6.8	
RDW-CV %	16.3	2.6		15.2	2.6		14.3	2.6	
MPV fL	10.3	1.6		9.9	1.1		9.8	0.9	

Refer to package insert for Sysmex recommendations on setting QC targets and limits
Distributed by
Sysmex America, Inc.