



**EIGHTCHECK-3WP X-TRA™**  
HEMATOLOGY CONTROL

poch-100i

Parameters	Low Abnormal			Normal			High Abnormal		
	Quality Control Data Due on the 20 <sup>th</sup> of each month			Quality Control Data Due on the 20 <sup>th</sup> of each month			Quality Control Data Due on the 20 <sup>th</sup> of each month		
	Expiration Date: 2010-05-19			Expiration Date: 2010-05-19			Expiration Date: 2010-05-19		
	Lot Number: 00400710			Lot Number: 00400711			Lot Number: 00400712		
	Mean	Limit	BCR	Mean	Limit	BCR	Mean	Limit	BCR
WBC x103/mL	2.7	0.4		6.5	0.7		16.9	1.7	
RBC x106/mL	2.37	0.21		4.46	0.27		5.41	0.27	
HGB g/dL	6.3	0.3		13.0	0.7		16.5	0.7	
HCT %	19.5	1.7		38.9	2.7		48.7	3.9	
MCV fL	82.3	13.1		87.2	10.7		90.0	11.1	
MCH pg	26.6	3.4		29.1	3.0		30.5	2.6	
MCHC g/dL	32.3	4.0		33.4	3.7		33.9	3.8	
PLT x103/mL	50	17		193	35		461	60	
Lymph (W-SCR) %	23.9	3.6		31.5	3.2		36.0	3.6	
Mixed (W-MCR) %	13.5	6.1		10.9	2.7		16.0	2.7	
Neut (W-LCR) %	62.6	5.0		57.6	3.5		48.0	3.8	
Lymph (W-SCC) x103/mL	0.6	0.2		2.0	0.5		6.1	0.7	
Mixed (W-MCC) x103/mL	0.4	0.2		0.7	0.2		2.7	0.5	
Neut (W-LCC) x103/mL	1.7	0.3		3.8	1.6		8.1	1.0	
W-SMV fL	67.3	16.8		66.1	14.5		66.0	13.2	
W-LMV fL	222.9	22.3		209.4	20.9		209.7	21.0	
RDW-SD fL	40.1	6.8		42.6	6.8		44.3	6.6	
RDW-CV %	13.8	2.2		13.9	2.4		13.9	2.5	
MPV fL	10.4	1.7		10.2	1.1		10.0	0.9	

Refer to package insert for Sysmex recommendations on setting QC targets and limits  
Distributed by  
Sysmex America, Inc.