



**EIGHTCHECK-3WP X-TRA™**  
HEMATOLOGY CONTROL

poch-100i

Low Abnormal				Normal				High Abnormal			
Data Due Date #1:		30-Aug-10		Data Due Date #2:		02-Oct-10		Data Due Date #3:		03-Nov-10	
Expiration Date:		2010-11-03		Expiration Date:		2010-11-03		Expiration Date:		2010-11-03	
Lot Number:		02080710		Lot Number:		02080711		Lot Number:		02080712	
Parameters	Mean	Limit	BCR	Mean	Limit	BCR	Mean	Limit	BCR		
WBC x10 <sup>3</sup> /μL	2.8	0.4		7.1	0.7		18.2	1.8			
RBC x10 <sup>6</sup> /μL	2.37	0.21		4.42	0.27		5.34	0.27			
HGB g/dL	6.3	0.3		13.0	0.7		17.0	0.7			
HCT %	19.3	1.6		38.3	2.7		49.1	3.9			
MCV fL	81.4	12.9		86.7	10.7		91.9	11.3			
MCH pg	26.6	3.4		29.4	3.0		31.8	2.7			
MCHC g/dL	32.6	4.0		33.9	3.8		34.6	3.8			
PLT x10 <sup>3</sup> /μL	47	16		203	37		470	61			
Lymph (W-SCR) %	24.8	3.7		31.2	3.1		35.8	3.6			
Mixed (W-MCR) %	9.2	4.1		10.5	2.6		15.3	2.6			
Neut (W-LCR) %	66.0	5.3		58.3	3.5		48.9	3.9			
Lymph (W-SCC) x10 <sup>3</sup> /μL	0.7	0.2		2.2	0.6		6.5	0.8			
Mixed (W-MCC) x10 <sup>3</sup> /μL	0.3	0.2		0.7	0.2		2.8	0.5			
Neut (W-LCC) x10 <sup>3</sup> /μL	1.8	0.4		4.2	0.6		8.9	1.1			
W-SMV fL	64.2	16.1		64.0	14.1		63.8	12.8			
W-LMV fL	204.9	20.5		204.7	20.5		204.7	20.5			
RDW-SD fL	40.8	6.9		41.9	6.7		45.0	6.8			
RDW-CV %	14.2	2.3		13.7	2.3		13.9	2.5			
MPV fL	10.1	1.6		9.9	1.1		9.7	0.9			

Refer to package insert for Sysmex recommendations on setting QC targets and limits  
Distributed by  
Sysmex America, Inc.