



EIGHTCHECK-3WP X-TRA™
HEMATOLOGY CONTROL

poch-100i

Parameters	Low Abnormal			Normal			High Abnormal		
	Mean	Limit	BCR	Mean	Limit	BCR	Mean	Limit	BCR
	Insight Due Date #1 8/1/2011			Insight Due Date #2 9/3/2011			Insight Due Date #3 10/5/2011		
	Expiration Date: 2011-10-05			Expiration Date: 2011-10-05			Expiration Date: 2011-10-05		
	Lot Number: 11790710			Lot Number: 11790711			Lot Number: 11790712		
WBC x10 ³ /μL	2.9	0.4		6.8	0.7		18.1	1.8	
RBC x10 ⁶ /μL	2.35	0.21		4.45	0.26		5.28	0.26	
HGB g/dL	6.0	0.3		12.8	0.7		16.6	0.7	
HCT %	18.2	1.6		38.0	2.7		47.7	3.9	
MCV fL	77.4	12.7		85.4	10.6		90.3	11.3	
MCH pg	25.5	3.4		28.8	3.1		31.4	2.7	
MCHC g/dL	33.0	4.1		33.7	3.8		34.8	3.9	
PLT x10 ³ /μL	49	17		194	36		445	60	
Lymph (W-SCR) %	25.0	3.7		29.7	3.1		37.3	3.5	
Mixed (W-MCR) %	10.7	4.6		13.1	2.8		15.9	2.7	
Neut (W-LCR) %	64.3	5.2		57.2	3.5		46.8	3.9	
Lymph (W-SCC) x10 ³ /μL	0.7	0.2		2.0	0.5		6.8	0.7	
Mixed (W-MCC) x10 ³ /μL	0.3	0.2		0.9	0.2		2.9	0.5	
Neut (W-LCC) x10 ³ /μL	1.9	0.4		3.9	0.9		8.4	1.1	
W-SMV fL	68.0	16.6		67.8	14.5		65.1	13.1	
W-LMV fL	201.8	20.5		203.4	20.3		204.3	20.3	
RDW-SD fL	41.5	7.1		42.8	6.8		45.0	6.8	
RDW-CV %	15.1	2.4		14.2	2.3		14.0	2.5	
MPV fL	10.0	1.6		9.6	1.1		9.5	0.9	

Refer to package insert for Sysmex recommendations on setting QC targets and limits
Distributed by
Sysmex America, Inc.