



EIGHTCHECK-3WP X-TRA™
HEMATOLOGY CONTROL

poch-100i

| Parameters | Low Abnormal | | | Normal | | | High Abnormal | | |
|-----------------------|-----------------------------|-------|-----|-----------------------------|-------|-----|-----------------------------|-------|-----|
| | Mean | Limit | BCR | Mean | Limit | BCR | Mean | Limit | BCR |
| | Data Due Date #1: 12-Aug-08 | | | Data Due Date #2: 09-Sep-08 | | | Data Due Date #3: 13-Oct-08 | | |
| | Expiration Date: 2008-10-08 | | | Expiration Date: 2008-10-08 | | | Expiration Date: 2008-10-08 | | |
| | Lot Number: 81830710 | | | Lot Number: 81830711 | | | Lot Number: 81830712 | | |
| | | | | | | | | | |
| WBC x103/mL | 2.5 | 0.4 | | 7.1 | 0.7 | | 18.1 | 1.8 | |
| RBC x106/mL | 2.38 | 0.21 | | 4.29 | 0.26 | | 5.19 | 0.26 | |
| HGB g/dL | 6.3 | 0.3 | | 12.6 | 0.6 | | 15.9 | 0.6 | |
| HCT % | 18.5 | 1.6 | | 35.9 | 2.5 | | 44.7 | 3.6 | |
| MCV fL | 77.7 | 12.3 | | 83.7 | 10.3 | | 86.1 | 10.6 | |
| MCH pg | 26.5 | 3.4 | | 29.4 | 3.1 | | 30.6 | 2.6 | |
| MCHC g/dL | 34.1 | 4.3 | | 35.1 | 3.9 | | 35.6 | 4.0 | |
| PLT x103/mL | 49 | 16 | | 187 | 34 | | 434 | 56 | |
| Lymph (W-SCR) % | 25.0 | 3.8 | | 29.7 | 3.0 | | 34.8 | 3.5 | |
| Mixed (W-MCR) % | 7.9 | 3.6 | | 12.2 | 3.1 | | 15.6 | 2.7 | |
| Neut (W-LCR) % | 67.1 | 5.4 | | 58.1 | 3.5 | | 49.6 | 4.0 | |
| Lymph (W-SCC) x103/mL | 0.6 | 0.2 | | 2.1 | 0.5 | | 6.3 | 0.8 | |
| Mixed (W-MCC) x103/mL | 0.2 | 0.1 | | 0.9 | 0.2 | | 2.8 | 0.5 | |
| Neut (W-LCC) x103/mL | 1.7 | 0.3 | | 4.1 | 0.6 | | 9.0 | 1.1 | |
| W-SMV fL | 65.3 | 16.3 | | 65.5 | 14.4 | | 65.1 | 13.0 | |
| W-LMV fL | 199.8 | 20.0 | | 193.4 | 19.3 | | 193.4 | 19.3 | |
| RDW-SD fL | 42.6 | 7.2 | | 42.7 | 6.8 | | 44.1 | 6.6 | |
| RDW-CV % | 15.5 | 2.5 | | 14.4 | 2.4 | | 14.5 | 2.6 | |
| MPV fL | 10.0 | 1.6 | | 9.8 | 1.1 | | 9.7 | 0.9 | |